



C.F.P.[®] Hip Prosthesis System



Presented by:



CE 0482

Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg, Germany

P.O. Box 63 05 52 · 22315 Hamburg, Germany

Tel.: +49 40 53995-0 · Fax: +49 40 5386929

E-mail: info@linkhh.de · www.linkorthopaedics.com

C.F.P.[®] Hip Prosthesis System

	System Description	
02	Bioharmony in Total Hip Replacement	
02	C.F.P. [®] Prosthesis Stem	
03	X-LINKed [®]	
03	Indications/Contraindications	
	Implants	
04	Prosthesis Stems	
05	Exchangeable Collars	
06	Prosthesis Heads A: BIOLOX [®] - Ceramic	
07	Prosthesis Heads B: CoCrMo	
	Instruments	
08	Instrument Set for LINK [®] C.F.P. [®] Prosthesis Stems	
13	Additional Instruments	
	Surgical Technique	
14	Preoperative Planning	
15	Surgical Approaches	
	Accessories	
24	X-ray Templates	
25	Literature	
27	Index	
	Important Information	

■ Bioharmony in Total Hip Replacement

The LINK® **C.F.P.®* Prosthesis Stem** and the **T.O.P.** Cementless Acetabular Cup System** have been designed especially for biologically young and active patients. With conventional hip prostheses such patients are subjected to increased rates of aseptic loosening because of their long life-expectancy.

The C.F.P.® and T.O.P. take into account the bio-mechanical load and fixation principles dictated by hip anatomy and physiology patterns, the retention of the femoral neck and uncemented fixation techniques to provide strong and stable implant anchorage. They thus provide an optimal basis for any future interventions.

* Collum Femoris Preserving

** Trabekel Orientierte Pfanne

■ C.F.P.® Hip Prosthesis Stem

- Minimal bone resection because femoral neck and proximal cancellous bone are preserved
- Parts of the circumflex artery's vascular network are preserved thereby maintaining the vascular supply of the femoral neck after resection of the femoral head
- Stem shape and surgical technique allow an uncemented implantation with up to 87% prosthesis-to-bone contact
- Prosthesis stem takes account of anatomical stem shape and anteversion
- Different stem curvatures enable secure support of the stem at the medial cortex
- Collar allows reintroduction of physiological loads into the femur
- Made of **Ti/otan**® (Ti6Al4V)

■ System Description – Indications/Contraindications

■ X-LINKed® – highly crosslinked Polyethylene for Acetabular Cups and Inserts

The LINK® Acetabular Cup Portfolio covers a wide range of indications for acetabular treatment. LINK® Acetabular Cups are used in many cases for both biologically young, active patients and older patients and are also used in surgical revision.

Prosthesis heads with a diameter of up to 36 mm can be used with cups made from X-LINKed® Polyethylene. This increases the range of movement and gives optimal joint stability.

By using X-LINKed® Polyethylene, abrasive wear can be reduced significantly which extends the service life of the implant components.

X-LINKed® Polyethylene is available for the following acetabular cup systems:

Cementless fixation

- T.O.P. Polyethylene inserts
- BetaCup® Polyethylene inserts

Cementable fixation

- Lubinus Eccentric/Dysplastic Polyethylene Acetabular Cup, without snap fit
- IP Polyethylene Acetabular Cup
- FAL Polyethylene Acetabular Cup, anti-luxation

■ Indications/Contraindications

Note: For specific indications/contraindications, see page 26.

■ Prosthesis Stems, cementless, taper 12/14 mm

Material: *Ti/zotan*® with HX® Coating (calcium phosphate)


HX® coated



Curvature A (slight)				Curvature B (strong)	
CCD angle 117°	CCD angle 126°	HX® coated		CCD angle 117°	CCD angle 126°
Item no.	Item no.	Size	Version	Item no.	Item no.
290-118/17	290-118/26	x-small	right	290-170/17	290-170/26
290-120/17	290-120/26	small	right	290-174/17	290-174/26
290-126/17	290-126/26	medium	right	290-178/17	290-178/26
	290-128/26	medium large	right		290-180/26
290-132/17	290-132/26	large	right	290-182/17	290-182/26
290-138/17	290-138/26	x-large	right	290-186/17	290-186/26
290-119/17	290-119/26	x-small	left	290-171/17	290-171/26
290-121/17	290-121/26	small	left	290-175/17	290-175/26
290-127/17	290-127/26	medium	left	290-179/17	290-179/26
	290-129/26	medium large	left		290-181/26
290-133/17	290-133/26	large	left	290-183/17	290-183/26
290-139/17	290-139/26	x-large	left	290-187/17	290-187/26

Material: *Ti/zotan*®

mikroporous ≈ 70 µm roughness



Curvature A (slight)				Curvature B (strong)	
CCD angle 117°	CCD angle 126°	mikroporous with 70 µm roughness		CCD angle 117°	CCD angle 126°
Item no.	Item no.	Size	Version	Item no.	Item no.
290-094/17	290-094/26	x-small	right	290-150/17	290-150/26
290-100/17	290-100/26	small	right	290-154/17	290-154/26
290-106/17	290-106/26	medium	right	290-158/17	290-158/26
	290-108/26	medium large	right		290-160/26
290-112/17	290-112/26	large	right	290-162/17	290-162/26
290-114/17	290-114/26	x-large	right	290-166/17	290-166/26
290-095/17	290-095/26	x-small	left	290-151/17	290-151/26
290-101/17	290-101/26	small	left	290-155/17	290-155/26
290-107/17	290-107/26	medium	left	290-159/17	290-159/26
	290-109/26	medium large	left		290-161/26
290-113/17	290-113/26	large	left	290-163/17	290-163/26
290-115/17	290-115/26	x-large	left	290-167/17	290-167/26

Curvature A or B			
Size	B mm	C mm	E mm
x-small	19.5	82	11
small	22	105	10
medium	24	115	11.5
medium large	25	120	12
large	26	125	13.5
x-large	28	135	15.5

■ Implants

■ Exchangeable Collars, cementless




All C.F.P.[®] hip prostheses are supplied with a standard size collar.

Material: Unalloyed titanium and HX[®] Coating (calcium phosphate)

Item no.	HX [®] coated, for stem size	Item no.
Standard		Large
290-190/01	x-small	290-191/01
290-190/02	small	290-191/02
290-190/03	medium	290-191/03
290-190/13	medium large	290-191/13
290-190/04	large	290-191/04
290-190/05	x-large	290-191/05


Material: Unalloyed titanium

Item no.	microporous surface, 70 µm roughness, for stem size	Item no.
Standard		Large
290-192/01	x-small	290-193/01
290-192/02	small	290-193/02
290-192/03	medium	290-193/03
290-192/13	medium large	290-193/13
290-192/04	large	290-193/04
290-192/05	x-large	290-193/05



Standard		For stem size	Large	
H mm	K mm		H mm	K mm
32	26	x-small	36	29
34	28	small	38	31
36	30	medium	40	33
37	31	medium large	41	34
38	32	large	42	35
40	34	x-large	44	37

HX[®] Coating (calcium phosphate)



The HX[®] Coating is an osteoconductive coating approximately 15 µm thick. The unique electrochemical coating process gives it immense mechanical strength, which is necessary for it to endure the stresses during implantation. The large porous structure provides the optimal environment for bone ongrowth.

The HX[®] Coating preserves the open-pored structure of the implant surface because it is so thin – thanks to the process of electrochemically coating.

■ **Prosthesis Heads**

Prosthesis heads A

Material: BIOLOX® forte* – ceramic



Item no.	Head Ø mm	Taper mm	Neck length mm
128-928/01	28	12/14	short -3.5
128-928/02	28	12/14	medium 0
128-928/03	28	12/14	long +3.5
128-932/01	32	12/14	short -4
128-932/02	32	12/14	medium 0
128-932/03	32	12/14	long +4
128-936/01	36	12/14	short -4
128-936/02	36	12/14	medium 0
128-936/03	36	12/14	long +4
128-940/01**	40	12/14	short -4
128-940/02**	40	12/14	medium 0
128-940/03**	40	12/14	long +4

Prosthesis heads A

Material: BIOLOX® delta* – ceramic



Item no.	Head Ø mm	Taper mm	Neck length mm
128-791/01	28	12/14	short -3.5
128-791/02	28	12/14	medium 0
128-791/03	28	12/14	long +3.5
128-792/01	32	12/14	short -4
128-792/02	32	12/14	medium 0
128-792/03	32	12/14	long +4
128-792/04**	32	12/14	extra long +7
128-793/01	36	12/14	short -4
128-793/02	36	12/14	medium 0
128-793/03	36	12/14	long +4
128-793/04**	36	12/14	extra long +8
128-794/01	40	12/14	short -4
128-794/02	40	12/14	medium 0
128-794/03	40	12/14	long +4
128-794/04**	40	12/14	extra long +8

All BIOLOX® forte* and BIOLOX® delta* components are inter compatible.

**On request

*BIOLOX® delta and BIOLOX® forte are made by CeramTec GmbH, Plochingen, Germany

■ Implants

Prosthesis heads B

Material: CoCrMo alloy



Item no.	Head Ø mm	Taper mm	Neck length mm
128-822/01	22	12/14	short -3.5
128-822/02	22	12/14	medium 0
128-824/01	24	12/14	short -3.5
128-824/02	24	12/14	medium 0
128-826/01	26	12/14	short -3.5
128-826/02	26	12/14	medium 0
128-826/03	26	12/14	long +3.5
128-828/01	28	12/14	short -3.5
128-828/02	28	12/14	medium 0
128-828/03	28	12/14	long +3.5
128-828/04**	28	12/14	extra long +10.5
128-832/01	32	12/14	short -4
128-832/02	32	12/14	medium 0
128-832/03	32	12/14	long +4
128-832/04**	32	12/14	extra long +8.5
128-836/01	36	12/14	short -4
128-836/02	36	12/14	medium 0
128-836/03	36	12/14	long +4
128-836/04**	36	12/14	extra long +8

**On request


■ Instrument Set for LINK® C.F.P.® Prosthesis Stems



Item no.	Instrument set, complete (container 1 + 2)
291-020/16	Set complete in 2 standard containers N11 & N31, on 3 trays with product illustrations and storage inserts
05-2003/03	N31 standard container, empty , 575 x 275 x 170 mm
05-2001/03	N11 standard container, empty , 575 x 275 x 100 mm
291-021/07	Upper tray (container 1), empty, perforated stainless steel, 550 x 265 x 50 mm
291-023/06	Lower Tray (container 1), empty, perforated stainless steel, 550 x 265 x 50 mm
291-022/07	Tray (container 2), empty, perforated stainless steel, 550 x 265 x 50 mm

■ Instruments


291-021/07 Upper Tray



Qty.

1	291-021/07	Upper tray (container 1), empty, perforated stainless steel, 550 x 265 x 50 mm	1
2	291-080	Guide to open femoral cavity, 180 mm	1
3	130-611	Impactor, 280 mm	1
4	291-081	Trocar awl, 200 mm	1
5	64-1181/06	Hex screwdriver, 175 mm, hex 2 mm	1
6	291-083/02	Stem inserter, 250 mm, 2 parts	1
7	130-393/60	Handle with quick coupling, 330 mm	1
8	130-609	Hex screwdriver, 175 mm, hex 4 mm	1
9	291-085	Stem extractor, 435 mm	1
10	291-082	Bone curette curved, to open femoral cavity, 260 mm	1
11	130-406/01A*	Calcar reamer Ø 40 mm, fitting for Jacobs Chuck, fitting optional*	1
	10-5371**	Hex screwdriver, 180 mm, hex 1.5 mm (**not included in instrument set)	1


*How to order: 131-406/01B = with Hudson fitting



A
Jacobs Chuck



B
Hudson



C
Harris



D
AO



H
Zimmer

System
Description

Indications/
Contraindications

Implants

Instruments


Surgical Technique

Accessories

Literature

Index

291-023/06 Lower Tray



		Qty.
12	291-023/06 Lower tray (container 1), empty, perforated stainless steel, 550 x 265 x 50 mm	1

Grey plastic trial heads with snap-fit for compressor cylinder

	Item no.	Ø	Neck length	Neck length mm	Qty.
13	291-826/01	26	short	-3.5	1
14	291-826/02	26	medium	0	1
15	291-826/03	26	long	+3.5	1
16	291-828/01	28	short	-3.5	1
17	291-828/02	28	medium	0	1
18	291-828/03	28	long	+3.5	1
	291-828/04*	28	extra long	+10.5	1
19	291-832/01	32	short	-4	1
20	291-832/02	32	medium	0	1
21	291-832/03	32	long	+4	1
	291-832/04*	32	extra long	+8.5	1
	291-832/05*	32	extra long	+7	1
22	291-836/01	36	short	-4	1
23	291-836/02	36	medium	0	1
24	291-836/03	36	long	+4	1
	291-836/04*	36	extra long	+8	1
	291-840/01*	40	short	-4	1
	291-840/02*	40	medium	0	1
	291-840/03*	40	long	+4	1
	291-840/04*	40	extra long	+8	1

*On request

■ Instruments

Curved probes with sizing olives, graduated, 300 mm

Item no.	Size	For prosthesis stem	Qty.
25 291-102/01	1	x-small	1
26 291-100/02	2	small	1
27 291-100/03	3	medium	1
28 291-100/35	35	medium large	1
29 291-100/04	4	large	1
30 291-100/05	5	x-large	1

Colored plastic trial heads, taper 12/14 mm



Item no.	Ø	Neck length	Neck length mm	Color	Qty.
31 131-926/01	26	short	-3.5	green	1
32 131-926/02	26	medium	0	blue	1
33 131-926/03	26	long	+3.5	black	1
34 131-928/01	28	short	-3.5	green	1
35 131-928/02	28	medium	0	blue	1
36 131-928/03	28	long	+3.5	black	1
131-928/04*	28	extra long	+10.5	brown	1
37 131-932/01	32	short	-4	green	1
38 131-932/02	32	medium	0	blue	1
39 131-932/03	32	long	+4	black	1
131-932/04*	32	extra long	+8.5	brown	1
131-932/05*	32	extra long	+7	orange	1
40 131-936/01	36	short	-4	green	1
41 131-936/02	36	medium	0	blue	1
42 131-936/03	36	long	+4	black	1
131-936/04*	36	extra long	+8	brown	1
131-940/01*	40	short	-4	green	1
131-940/02*	40	medium	0	blue	1
131-940/03*	40	long	+4	black	1
131-940/04*	40	extra long	+8	brown	1


Trial neck sections

Item no.	CCD angle	Version	Qty.
43 291-120/17	117°	right	1
44 291-121/17	117°	left	1
45 291-120/26	126°	right	1
46 291-121/26	126°	left	1

47 130-600	Driver for prosthesis heads with exchangeable plastic head, 170 mm	Qty. 1
----------------------------	--	-----------

*On request

291-022/07 Tray



48 291-022/07 Tray (container 2), empty, perforated stainless steel, 550 x 265 x 50 mm 1 Qty.

Bone compressors, curvature A

Item no.	Size	Version	Item no.	Size	Version
49 291-200/03	x-small	right	55 291-200/04	x-small	left
50 291-201/01	small	right	56 291-201/02	small	left
51 291-202/01	medium	right	57 291-202/02	medium	left
52 291-202/11	medium large	right	58 291-202/12	medium large	left
53 291-203/01	large	right	59 291-203/02	large	left
54 291-204/01	x-large	right	60 291-204/02	x-large	left

Bone compressors, curvature B

Item no.	Size	Version	Item no.	Size	Version
61 291-210/03	x-small	right	67 291-210/04	x-small	left
62 291-211/01	small	right	68 291-211/02	small	left
63 291-212/01	medium	right	69 291-212/02	medium	left
64 291-212/11	medium large	right	70 291-212/12	medium large	left
65 291-213/01	large	right	71 291-213/02	large	left
66 291-214/01	x-large	right	72 291-214/02	x-large	left

■ Instruments

■ Additional Instruments

(not included in instrument set for C.F.P.® Prosthesis Stems)



130-601

Replacement head for driver 130-600



130-165

Mallet

Ø 30 mm, 270 mm, 600 gram

130-393/15

**Handle for rasp stems,
right hip, angled**



130-393/25

**Handle for rasp stems,
left hip, angled**

■ Preoperative Planning

The aim in preoperative planning is to establish the approximate size of implant required and the optimal position in which to place it. Careful planning helps surgeons to foresee and avoid surgical problems.

For the best possible results the appropriate implant should be selected using C.F.P.[®] X-ray templates which are available at a scale of 1.1:1. When used in combination with recent pelvic X-rays (A/P and M/L views) these templates serve as a useful aid in planning surgical procedure and determining implant size.


In planning the resection level, the femoral neck must be considered along with centre of rotation and leg length. The femoral neck should remain intact as far as possible so that the original anatomy can be reconstructed.

The choice of implant should ensure that the greatest possible length of hip stem lies along the Adam's Bow. Stems are available with two different curvatures, curve A and curve B, for this purpose.


The distal tip of the prosthesis should run along the centre of the medullary canal and should not have any distal cortical contact.

■ Surgical Technique


■ Surgical Approaches



A: **Watson Jones**



B: **Hardinge**





C: **Moore**

The choice depends on the surgeon's experience and his/her decision based on the individual situation.

The following approaches are common:

- antero-lateral - **Watson Jones** (A)
- lateral - **Hardinge** (B)
- postero-lateral - **Moore** (C)





1



Dislocation of Femoral Head

After dislocation of the femoral head the femoral neck and the proximal rim of both trochanters are exposed and existing osteophytes of the femoral head removed.

2 and 3

Femoral Head Resection

During preoperative planning the osteotomy site is determined by applying tangents as shown in Fig. 3. The line of the resection on the isthmus corresponds to the line connecting the intersection points of the tangents. The distance from the base of the major trochanter is usually 1.5 cm.







■ Surgical Technique

3

Preparation of Proximal Femur

Using a trocar awl a small hole is made in the resected femoral neck, as far medially as possible, to accommodate the pin of the guide.



4

Targeting opening of femoral canal


The guide is adjusted to the planned prosthesis size (A) and placed on the resection site taking the premarked hole into consideration. The trocar awl is used to site the opening of the medullary cavity.



5

Opening of Femoral Canal

The femoral canal is opened using a curved bone curette.




6

Determination of Stem Size

To determine the stem size a curved probe is introduced into the canal starting with one size below the size planned. Each curved probe is inserted up to the ring mark.

The section below to the ring mark corresponds to the length of the bone compressor and of the final implanted stem.





7

Placing of Compressor


Finally the femoral canal is prepared for seating of the prosthesis stem by compressing the cancellous bone with the bone compressor. The choice of bone compressor depends on the stem curvature determined at the planning stage and the size indicated by the curved probe.

Initially a compressor should be used which is one or two sizes smaller than the planned size of the final implant. At this stage it is important to ensure that the medial teeth of the compressor do not remove the cortical bone at the calcar bow. When the smallest C.F.P.[®] prosthesis stem is used preparation of the femoral canal with a smaller compressor is not necessary. In this case preparation of the femoral bone with a curved bone curette and the sizing olive is sufficient.

8


Fitting of Compressor

The bone compressor is driven in below the resection level so that reaming with the calcar reamer produces a smooth surface on which the removable collar can be seated.



Attention:

To prevent the reamer from being damaged it must always be pushed as far as possible onto the guide pin before starting to ream.



9


Preparation of Trial Reduction

Trial neck (left or right and 117° or 126° CCD angle) and trial head are attached to the bone compressor in situ. The bone compressor now serves as a trial stem.

10

Trial Reduction

The trial reduction indicates the correct position and the head-neck length.






11

Implantation of Stem

The final C.F.P.® stem is placed in the stem inserter. To prevent damage to the taper of the stem it is important to ensure that the inserter locks firmly onto the end face of the taper and holds the stem correctly beneath the taper. The C.F.P.® stem is then inserted as far as possible into the femoral canal using the stem inserter.

The design of the C.F.P.® stem incorporates the anatomical anteversion of the femoral neck. Adjustment of the anteversion, as performed when implanting straight stems, would lead to rotational malpositioning and must therefore be avoided.

Anteversion




12

Driving in Stem

Placing the tip of the impactor into the recess on the lateral collar the C.F.P.® stem can be driven home into the canal. Before it is finally driven home the collar may be removed in order to insert additional bone into the grooves.

The fixation screw on the collar must then be firmly tightened anticlockwise.



Note: LINK strongly recommends to not remove the collar. In case the surgeon does not follow this advice and the stem is implanted without collar, he/she has to make sure that the collar locking screw is removed as well.

13

Final Trial Reduction with Trial Head


As a precaution, a final trial run is performed using colored plastic trial heads.




14

Positioning of Prosthesis Head

The final prosthesis head is then mounted on the carefully cleaned taper of the stem and fixed with a light blow on the impactor.



■ Surgical Technique



15 and 16

The C.F.P.® Prosthesis Stem in situ.
Final trial reduction with permanent
implant components.

System
Description

Indications/
Contraindications

Implants

Instruments

Surgical Technique

Accessories

Literature

Index

■ X-ray Templates

X-ray templates for C.F.P.[®] Prosthesis Stems, cementless (with neutral head-neck length), taper 12/14 mm, 110% actual size

Curvature A			
Item no.	Head Ø mm	CCD angle	Stem size
290-258/52	28, 32, 36, 40	117°	x-small, small, medium, large, x-large
290-258/32	28, 32, 36, 40	126°	x-small, small, medium, medium large, large, x-large
Curvature B			
Item no.	Head Ø mm	CCD angle	Stem size
290-259/52	28, 32, 36, 40	117°	x-small, small, medium, large, x-large
290-259/32	28, 32, 36, 40	126°	x-small, small, medium, medium large, large, x-large

Instructions for Cleaning and Maintenance

Specific instructions for instruments are available on request via E-mail customer@linkhh.de

Literature

F. Pipino, L. Molfetta, M. Grandizio

Preservation of the femoral neck in hip arthroplasty: results of a 13- to 17-year follow-up
Orthopaedics and Traumatology 2000 1:31-39 **(H110)**

F. Pipino

CFP prosthetic stem in mini-invasive total hip arthroplasty
Orthopaedics and Traumatology 2004 4:165-171 **(H137)**

F. Biggi, F. Pipino, F. Franclin, R. Lovato

DEXA evaluation of total hip arthroplasty with neck-preserving technique: 4-year follow-up
Orthopaedics and Traumatology 2004 5:156-159 **(H138)**

1. F. Pipino

1. Tissue sparing surgery (TSS) in hip and knee arthroplasty
 2. Tissue-sparing surgery: 25 years' experience with femoral neck preserving hip arthroplasty
Orthopaedics and Traumatology 2006 Vol. 7, Nr. 1 **(H144)**

2. F. Pipino, A. Keller

S. M. Röhrli, M. G. Li, E. Pedersen, G. Ullmark, B. Nivbrant

Migration Pattern of a short Femoral Neck Preserving Stem
Clinical Orthopaedics & Related Research 2006 Nr. 448, pp. 73-78 **(H146)**

H. Wacha, G. Domsel, R. Motz

Follow-up-Studie – 3 Jahre nach Implantation einer Collum femoris erhaltenden Totalendoprothese der Hüfte
Orthopädische Praxis Jahrg. 43, Heft 8 2007 **(H151)**

G. von Foerster

Minimal Invasive Hüft-Endoprothetik posteriorer Zugang – langzeitiger Nutzen?
MOT 5/2007 **(H153)**



Link News - Orthopedics Update

Case Histories:
 C.F.P.® Hip Prosthesis Stem
 T.O.P.® Acetabular Cup
 Order number: Nr.10 en (English)

Link News Spezial - Orthopädie aktuell

Veranstaltung Frankfurt a. M., 2002
 C.F.P.® Hüftprothesenschaft
 T.O.P.® Hüftpfanne
 Order number: Nr.15 dt (German)

Link News - Orthopädie aktuell

Knochensparende Wechseloperation mit dem C.F.P.®-Hüftprothesenschaft
 Order number: Nr.19 dt (German)


D. Briem, M. Schneider, N. Bogner, N. Botha, M. Gebauer, T. Gehrke, B. Schwantes

Mid-term results of 155 patients treated with a collum femoris preserving (CFP) short stem prosthesis
International Orthopaedics (SICOT), July 2010

M. Pons

Learning curve and short-term results with a short-stem C.F.P.® system
HIP International, Vol. 20, Suppl. 7, pp. 52-57, April-June 2010

Catalog
 T.O.P.
 Order number:
 672en



Catalog
 BetaCup®
 Order number:
 603en



Catalog
 MIT-H®
 Order number:
 674en

	Hip Prosthesis	Prosthesis Heads	
	C.F.P. [®] Total Hip Prosthesis	BILOX [®] forte* + delta* Prosthesis Heads A	Prosthesis Heads B
General Indications			
Mobility-limiting diseases, fractures or defects which cannot be treated by conservative or osteosynthetic procedures	X	X	X
Indications			
Primary and secondary coxarthrosis	X	X	X
Osteoarthritis	X	X	X
Necrosis of the femoral head	X	X	X
Femoral neck fractures		X	X
Revision after implant loosening		X	X
Contraindications			
Poor general state of health	X	X	X
Acute and chronic infections, local and systemic	X	X	X
Allergies due to (implant) materials	X	X	X
Distinctive muscular-, nerve-, vascular or other diseases which put the affected limb at risk	X	X	X
Insufficient/inadequate bone mass which prevents a stable anchor of the prosthesis	X	X	X
Relative Contraindications			
Adiposity	X	X	X
Lacking or foreseeable not assured compliance	X	X	X
Foreseeable overload/overstressing of the joint prosthesis	X	X	X
Osteoporosis	X		

Please note: These indications/contraindications refer to standard cases. The ultimate decision on whether or not an implant is suitable for a patient must be made by the surgeon based on his/her individual analysis and his/her experience.

*BILOX[®] forte and BILOX[®] delta are made by CeramTec GmbH, Plochingen, Germany

■ Index

	page		page
05-2001/03	08	290-094/17 to 290-187/17	04
05-2003/03	08	290-190/01 to 290-190/13	05
		290-191/01 to 290-191/13	05
10-5371	09	290-192/01 to 290-192/13	05
		290-193/01 to 290-193/13	05
64-1181/06	09	290-258/32, 290-258/52	24
		290-259/32, 290-259/52	24
128-791/01 to 128-791/03	06	291-020/16	08
128-792/01 to 128-792/04	06	291-021/07	08, 09
128-793/01 to 128-793/04	06	291-022/07	08, 12
128-794/01 to 128-794/04	06	291-023/06	08, 10
128-822/01, 128-822/02	07	291-080	09
128-824/01, 128-824/02	07	291-081	09
128-826/01 to 128-826/03	07	291-082	09
128-828/01 to 128-828/04	07	291-083/02	09
128-832/01 to 128-832/04	07	291-085	09
128-836/01 to 128-836/04	07	291-100/02 to 291-100/35	11
128-928/01 to 128-928/03	06	291-102/01	11
128-932/01 to 128-932/03	06	291-120/17, 291-120/26	11
128-936/01 to 128-936/03	06	291-121/17, 291-121/26	11
128-940/01 to 128-940/03	06	291-200/03, 291-200/04	12
		291-201/01 to 291-204/01	12
130-165	13	291-201/02 to 291-204/02	12
130-393/15, 130-393/25	13	291-210/03, 291-210/04	12
130-393/60	09	291-211/01 to 291-214/01	12
130-406/01A-H	09	291-211/02 to 291-214/02	12
130-600	11	291-826/01 to 291-826/03	10
130-601	13	291-828/01 to 291-828/04	10
130-609	09	291-832/01 to 291-832/05	10
130-611	09	291-836/01 to 291-836/04	10
		291-840/01 to 291-840/04	10
131-926/01 to 131-926/03	11		
131-928/01 to 131-928/04	11		
131-932/01 to 131-932/05	11		
131-936/01 to 131-936/04	11		
131-940/01 to 131-940/04	11		

Please note the following regarding the use of our implants:

1. Choosing the right implant is very important.

The size and shape of the human bone determine the size and shape of the implant and also limit the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.

2. Correct handling of the implant is very important.

Under no circumstances should the shape of a finished implant be altered, as this shortens its life span. Our implants must not be combined with implants from other manufacturers.

The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.

3. Implants must not be reused.

Implants are supplied sterile and are intended for single use only. Used implants must not be reused.

4. After-treatment is also very important.

The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of healthy bone!

5. Unless otherwise indicated, implants are supplied in sterile packaging.

Note the following conditions for storage of packaged implants:

- Avoid extreme or sudden changes in temperature.
- Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the "Use by" date indicated on the packaging.
- They must not be exposed to frost, dampness or direct sunlight, or mechanical damage.
- Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The "Use by" date is indicated on the product label.
- Do not use an implant if the packaging is damaged.

6. Traceability is important.

Please use the documentation stickers provided to ensure traceability.

7. Further information on the material composition is available on request from the manufacturer.

Follow the instructions for use!

Waldemar Link GmbH & Co. KG, Hamburg

All content in this catalog, including text, pictures and data, is protected by copyright. Every instance of use, whether in part or in whole and which is not permitted by the copyright act, is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, publishing, saving, processing, or passing on of content stored in databases or other electronic media and systems, in any manner or form. The information in the catalogs is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his/her responsibility to duly consider the particularities of each individual case.



Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg, Germany
P.O. Box 63 05 52 · 22315 Hamburg, Germany
Tel.: +49 40 53995-0 · Fax: +49 40 5386929
E-mail: info@linkhh.de · www.linkorthopaedics.com

LINK[®]

